

# Good Laboratory Work Practices (GLWP)

Nazlina Ibrahim & Clarence Ko







Sandia National Laboratories

## **Synopsis**

Good Laboratory Work Practices (GLWP) are practices and procedures that can reduce or mitigate biorisk. GLWP, along with Personal Protective Equipment (PPE) and Engineering Controls, is a part of risk mitigation







# Scope

- Good housekeeping
- Biosafety cabinet
- Proper pipetting skill
- Hand hygiene
- Field GLWP (optional)







### Learning Objectives

- Participants will be exposed to the importance of following GLWP
- Participants will understand some common good laboratory work practices
- Participants will recognize tasks within the laboratory's biosafety level with exposure to hazards








# Review the implementation of biorisk management system

- ► Good Laboratory Work Practices (GLWP) are applicable to mitigation in the AMP model
- Evaluate and ensure the system is working according to the risk management goal









# The AMP Model: Recognize suitable GLWPs according to the risk

- Plan-Do-Check-Act Cycle aligned with AMP Model for biorisk management
- Plan: Planning, identification of hazards and risk, and establishing goals
- ▶ Do: Implementing training and operational issues
- Check: Checking, monitoring, and corrective action
- Act: Reviewing, process innovation, and acting to make needed changes to the management system
- ► Taken together, the three elements of AMP constitute a complete biorisk management


The elements of the AMP model also underpinned

CWA 15793:2011 - Laboratory Biorisk Management Standard. (ISO 35001 - Laboratory Biorisk Management Standard, 2019)











# Key Components of Biorisk Management in GLWP

Identify hazard

Evaluate risk

Identify mitigation plans

Provide/ Follow guideline for GLWP










## Hierarchy of Controls

Where can we position GLWP in the hierarchy of controls?











### Good Laboratory Work Practices

- Practices and procedures used in the laboratory when followed protect laboratory workers and the environment from the risk of exposure or reduce the risk of hazardous agents
- Why do we practice GLWP? Following GLWP ensures:
  - own safety,
  - quality of experimental results, and
  - organization overall safety








#### GLWP is practiced to avoid:

Laboratory Acquired Infections (LAIs)













(Determined through biological risk assessments, evaluated by the study of infectious agents or toxins, their transmission, and their ability to cause infection)









### Common sense in Laboratory Safety

- Laboratory has many hazards e.g. source chemical, biological, electrical, ergonomic, fire
- Participants must be aware of these hazards and know how to handle them appropriately
- Participants must pay attention to these hazards to ensure a safe laboratory environment

https://www.creative-biolabs.com/blog/index.php/common-sense-in-laboratory-safety/









#### **ACTIVITY:**

# LIST RELEVANT GOOD LABORATORY PRACTICES THAT REQUIRE COMMON SENSE

Can you think of practices that require common sense?

Write one practice per sticky note You have 5 minutes to answer this question









#### Common Sense in Laboratory Safety

#### Some examples

- chemical:
  - when handling strong acids and alkalis, avoid spills
  - do not touch these chemicals with bare hands, put on/use appropriate equipment when handling chemicals

Can you think of other practices that require common sense?

#### fire/flame:

- alcohol or combustible liquids must be poured away from fire sources
- flame on culturing glass rod must first be extinguished before dipping in methanol







https://www.creative-biolabs.com/blog/index.php/common-sense-inlaboratory-safety/

### **Poor Laboratory Work Practice**

- Expired chemicals still in the laboratory
- Medium agar (with or without culture) left until dry in the fridge - not stocked in a proper manner
- Plate in biohazard plastic and not autoclaved
- Working alone during weekends
- Switching on the radio/audio equipment loudly

ANYMORE THAT YOU KNOW ???????







#### Some Common GLWP

- Good housekeeping:
  - Daily decontamination of work surfaces and use of chemicals with decontamination properties, e.g. 70% ethanol,0.5% sodium hypochlorite, or 70% isopropanol (depending on the type of pathogens and work surface)
  - Use of an autoclave to reduce contamination
  - use biosafety cabinet if needed



Source: https://www.mynewlab.com/blog/laboratory-equipment-maintenance-101/







#### Biological safety cabinet (BSC)

- Biological safety cabinets are designed for the safe handling of pathogens and other potentially hazardous materials with minimal risk
- Proper use of a biosafety cabinet reduces the laboratorian's risk of exposure to potential pathogens and minimizes contamination of the work environment







#### Some Common GLWP

- Biosafety cabinet
  - type of BSC (according to the pathogen being handled)
  - technical competency must be acquired to use BSC properly



Source: https://www.freepik.com/premium-photo/scientist-work-sterile-hood\_27067196.









# Good Laboratory Work Practice related to Biosafety Cabinets



Watch and list 5
Good Work Practices
related to Biosafety
that you can identify
in the video



https://www.youtube. 2com/watch?v=96aZLom340&t=11s









# Poor Laboratory Work Practice - Activity

Watch and list 5 Poor Work Practices related to Biosafety that you can identify in the video

(Duration 3.11 minutes)







https://www.youtube.com/watch?v=dAOTg-3ev-w







# Biosafety mistakes that can be identified in the video

- Hair not covered. Hair must be tied up to avoid the possibility of dipping into chemicals/solutions
- Touching hair with gloved and ungloved hands while handling biological materials
- Biosafety cabinet filled with a lot of stuff that reduces the airflow
- The BSC airflow obstructed by items
- Presence of a non-researcher in the laboratory
- Use of soiled gloves
- Use of improper gloves
- Not wearing shoes
- Not putting on the laboratory gown properly
- Reusing the microcentrifuge







#### Some Common GLWP

- Proper pipetting skill is important
  - To avoid the formation of aerosol
  - For volume accuracy
  - Use of proper equipment (no mouth pipetting)









Source: Photo by Thirdman from Pexels: https://www.pexels?com/photo/a-scientist-using-a-laboratory-equipment-8940510/



#### Potential Aerosol Generating Activities in Laboratory



Source: https://www.artel.co/learning\_center/best-practices-use-micropipets/



Source: https://www.medical-supply.ie/2019/04/bunsen-burner-biological-cabinets/



Source: https://blog.microbiologics.com/7-streak-plate-method-best-practices/







Source: h National ws-medical.net/Hand-Held-Cordless-Homogenizer-System



Source: https://www.sciencephoto.com/media/607032 /view/blood-sample



# Proper Pipetting Skill

- Avoid errors with wrong pipetting technique (volume error)
- Avoid the formation of aerosol that contains Microbes
- Pre-wet the pipette tip prior to initial delivery to avoid sample loss due to evaporation, especially when handling volatile solutions such as organic solvents.
- No tip wiping
- Choose the correct pipet tips
- Do not pipet at an angle
- No mouth pipetting- use a device












#### Aerosols in the Lab

- Aerosols are colloidal systems of droplets and/or solid particles suspended in the air
- Being small in size (< 10 μm), they can remain suspended in the air for several hours and get widely dispersed by the air currents. If inhaled, they can cause severe infections
- Aerosolization takes place through breathing, speech, coughing, and sneezing. Common laboratory activities such as centrifugation, pipetting, opening of ampoules, and shaking can also cause aerosols







http://www.eppendorf.com/stayinformed







### Prevention against Aerosols in the Lab

- 1. Use aerosol-tight caps/lids during centrifugation
- 2. Consider the min./max. filling volume of vessels
- 3. Do not exceed the centrifuge speed limit for the vessel
- 4. Use plastic vessels: glass vessels are more likely to break or leak.
- 5. Alternatively, open and load/unload the rotor in a biosafety cabinet
- 6. If a tube breaks or leaks, wait 30 minutes before opening the centrifuge or aerosol-tight lid or cap, to allow for the aerosol settle



## WHY IT IS IMPORTANT TO PREVENT AEROSOL FORMATION?

- Many documented cases of laboratoryacquired infections result from the inhalation of infectious aerosols
- ➤ These infections may be more or less severe, depending on the agents, and present a health hazard for everyone in the laboratory

(\*1979-2005: Harding, L.H. and K.B. Byers. 2006. Epidemology of laboratory-acquired infections.)









# Hand hygiene

Will be discussed more in PPE module







Sandia National Laboratories

#### Some Common GLWP

- Hand hygiene
  - before and after the use of the laboratory
  - use of gloves is part of the principle in achieving hand hygiene



Source: https://www.freepik.com/photos/washing-hands". Washing hands photo created by freepik - www.freepik.com











Source:

https://www.infosihat.gov.my/pener

bitan-

multimedia/poster.raw?task=callele ment&item\_id=7420&element=3f919

76b-5428-4aba-b4c9-

d843a3316027&method=download&args[0]=0











Source:

https://www.infosihat.gov.my/penerbitan-

multimedia/poster.raw?task=callele ment&item\_id=7420&element=3f919 76b-5428-4aba-b4c9-

d843a3316027£method=download&a rgs[0]=0









# Good Work Practices (GWP) in the Field







#### GWP in the Field

- Remember that hazards are everywhere
- Do your AMP when planning to go on a field trip
- ▶ The goal of GWP in the field are:
  - for participants to return from the trip safely
  - for productive and valuable research and learning experience
  - to return samples, related equipment, and supplies in Good Working Order









#### PRIMARY PROBLEMS IN THE FIELD

- ► Injury: Trip and Fall Hazards
- Improper clothing (poor footwear, etc)
- Heavy Work Overload: 12 hours of fieldwork in 8 hours
- Heat Stroke, Heat Exhaustion
- Hypothermia
- Dehydration
- Sun burn
- Contact with wildlife, leeches, ticks, snakes, bacteria, etc.
- Automobile Incidents
- The public
- Getting Lost: compass, local map, GPS, and road map
- Drowning: sea, river, flooding (undercurrents, water column, excessive river flow)







#### PREPARE FOR UNCERTAIN CONDITIONS

- Know the work conditions and the environment
- Know what specific safety equipment is needed?
- Know the work hazards
- Bring a field first-aid kit
- Inform at least two people about where you are going and when you will return
- Bring water, a flashlight, matches, and a snack Just in Case!
- Check the weather forecast to prepare the necessary work clothing for the condition
- Work in Pairs
- For all fieldwork, field clothing/PPE will include long panties and hiking boots

# Suggested GLWP when returning from field

- Decontamination/Quarantine of plant or animal samples: reduce/eliminate the possibility of bringing spores into the laboratory/building
- Take precautions when handling
  - animal feces that could contain microbes, which may be dispersed by aerosol
  - animal tissues/body fluid that may contain infectious microbes







#### Field etiquette

- Leave the area as you found it
- Respect private property: This is a good reason to have a camera handy to record evidence in case of issues
- Bring your trash home: Do not bury or leave your trash or waste in the field
- Fill in holes, test pits, auger holes, percolation holes, shovel pits, and cores
- When walking keep a distance of at least 6 feet between people







### Key Messages

- 1. Good Laboratory Work Practices are techniques and methods of doing work in the laboratory that reduces biorisk.
- 2. Practicing GLWP is to ensure the safety of people working in the laboratory and others in the vicinity
- 3. Barriers to Good Laboratory Work Practices can be overcome using various strategies
- 4. Practice GLWP and make it a habit







# Thank you

Any question?





